

## Levapren® product range

	Vinyl acetate content, weight (%)	Viscosity ML (1+4) at 100 °C
<b>Standard grades:</b>		
Levapren® 400	40 ± 1.5	20 ± 4
Levapren® 450	45 ± 1.5	20 ± 4
Levapren® 500 HV	50 ± 1.5	27 ± 4
Levapren® 600 HV	60 ± 1.5	27 ± 4
Levapren® 650 HV VP*	65 ± 1.5	27 ± 4
Levapren® 700 HV	70 ± 1.5	27 ± 4
Levapren® 800 HV	80 ± 2.0	28 ± 6
Levapren® 900 HV	90 ± 2.0	approx. 38
<b>Pre-crosslinked grades:</b>		
Levapren® VP 500 XL*	50 ± 1.5	approx. 55
Levapren® VP 600 XL*	60 ± 1.5	approx. 55
Levapren® VP 700 XL*	70 ± 1.5	approx. 60
Levapren® VP 800 XL*	80 ± 2.0	approx. 55
Levapren® VP KA 8784*	70 ± 1.5	approx. 60
Levapren® VP KA 8815*	60 ± 1.5	approx. 55

\* Trial Product VP = Versuchsprodukt

## Levapren® packaging

	Pallet, 40 sacks 25-kg PE sack in a box	Pallet, 40 sacks 25-kg PE sack	Pallet, 35 sacks 25-kg PE sack
<b>Standard grades:</b>			
Levapren® 400		●	
Levapren® 450		●	
Levapren® 500 HV		●	
Levapren® 600 HV	●		
Levapren® 650 HV VP*	●		
Levapren® 700 HV	●		
Levapren® 800 HV	●		
Levapren® 900 HV	●		
<b>Pre-crosslinked grades:</b>			
Levapren® VP 500 XL*	●		
Levapren® VP 600 XL*	●		
Levapren® VP 700 XL*	●		
Levapren® VP 800 XL*	●		
Levapren® VP KA 8784*			●
Levapren® VP KA 8815*			●

\* Trial Product VP = Versuchsprodukt